

Heating Mantle 1000 ml, 300 W

MODEL - BT-018/D

MAKE - AK LABS

This heating mantle offers sturdy construction, dependable heating performance, and a durable insulated body for regular laboratory use. Its well-finished design provides reliable support for compatible flasks, while the adjustable control and clear indicator ensure convenient, consistent, and practical operation.

PRODUCT SPECIFICATION

Capacity	1000 ml
Wattage	300 W
Control Range	0-100
Heating Type	Electric Heating Mantle
Flask Compatibility	Round-Bottom Flask
Indicator	Main Power Indicator Lamp
Chamber	Insulated Heating Chamber

PRODUCT DETAILS

This heating mantle is suitable for heating liquids in compatible 1000 ml round-bottom flasks during evaporation, dissolution, reaction, and sample preparation. Its 300 W heating element provides efficient heat output, while the adjustable 0-100 control enables flexible heating regulation. The insulated chamber offers practical protection, and the built-in indicator lamp provides clear visual confirmation of power operation during laboratory procedures.