

Melting Point Apparatus [Mechanical Type (Heavy Type) Die-Cast]

MODEL - BT-033/B

MAKE - AK LABS

The Mechanical Type Heavy Duty Melting Point Apparatus features a robust die-cast body engineered for durability and dependable laboratory performance. Its solid construction, quality finish, clearly marked controls, and sturdy design make it suitable for repeated use in demanding laboratory environments.

PRODUCT SPECIFICATION

Type	Mechanical, Heavy Duty
Construction	Die-Cast Body
Heating Control	Variable 0-100
Heater Setting	Low / High
Sample Observation	Viewing Aperture
Indicator	Load Indicator
Operation	Electric Heating

PRODUCT DETAILS

Designed for determining the melting point of solid substances, this apparatus provides controlled heating through an adjustable control dial and separate heater settings. The viewing aperture allows observation of the sample during heating, while the illuminated load indicator shows operating status. Its heavy die-cast construction offers stability and durability, making it suitable for educational, research, pharmaceutical, chemical, and quality-control laboratories.