

Test Tube Holder (Wooden Handle)



MODEL - CH-123

MAKE - AK LABS

AK Labs Test Tube Holder is a laboratory-grade holder featuring a durable wooden handle and heat-resistant gripping jaws. Available in Iron Chrome Plated and Brass construction, it ensures safe handling of hot test tubes during heating, chemical experiments, and laboratory demonstrations.

PRODUCT SPECIFICATION

Product Type	Test Tube Holder with Wooden Handle
Material Options	Iron Chrome Plated / Brass
Handle Material	Seasoned Wooden Handle
Grip Mechanism	Spring Action Clamp
Finish	Chrome Plated / Polished Brass
Construction	Superior Laboratory Quality
Applications	Heating Test Tubes, Chemical Experiments, Laboratory Demonstrations & Educational Use

PRODUCT DETAILS

Manufactured from premium-quality Iron Chrome Plated or Brass, the AK Labs Test Tube Holder provides a secure grip for handling heated test tubes safely. The ergonomically designed wooden handle offers excellent heat insulation and comfortable handling during experiments. Ideal for chemistry laboratories, educational institutions, research facilities, and routine laboratory applications requiring safe manipulation of glassware over direct flame.