

Screw Gauge (Micrometer), Stainless Steel Rod

MODEL - PH-300

MAKE - AK LABS

AK Labs Screw Gauge with Stainless Steel Rod is a precision measuring instrument designed for accurate measurement of small diameters and thicknesses. Manufactured with a durable stainless steel spindle and 0.01 mm least count, it delivers reliable performance for laboratory and educational applications.

PRODUCT SPECIFICATION	
Product Type	Screw Gauge (Micrometer)
Brand	AK LABS
Measuring Rod	Stainless Steel
Least Count	0.01 mm (1/100 mm)
Construction	Hardened Steel Frame with Stainless Steel Spindle
Measuring Range	20 mm / 25 mm
Applications	Physics Labs, Engineering Workshops, Schools, Colleges & Industrial Measurement

PRODUCT DETAILS
The AK Labs Screw Gauge (Micrometer) is fitted with a precision-ground stainless steel measuring rod for long-lasting accuracy and corrosion resistance. Featuring a least count of 0.01 mm (1/100 mm), it is ideal for measuring wire diameter, sheet thickness, and small engineering components. Suitable for physics laboratories, engineering workshops, schools, colleges, and industrial quality control.